

## SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)

Ref. 223004



Flexible and reliable protection to extend the lifespan of your Master Palm 2 or 3 blood pressure monitor.

### Description

Silicone protection for MASTER PALM 2 and MASTER PALM 3 oximeters.  
Set of 2 protections, one yellow and one transparent.

29/10/2025 01:58:12 – For information purposes only, can be subject to change

## SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)

Ref. 223004



### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | NONE          | Medical device            | No       |
| Bar code                    | 3700446025434 | HS Code                   | 90330000 |
| Country                     | China         |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.074    |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Packed unit / Quantity      | 1             |                           |          |