

PAIR OF PINK EARTIPS FOR PULSE II

Ref. 506904



Comfortable and colorful ear tips designed for standard use with the Pulse II.

Description

Pair of ear tips for SPENGLER PULSE II. Standard size, pink color.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030025	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		