

**NAVY BLUE CHROME PLATED BRASS TUBE
BINAURALS****Ref. 507912**

navy blue chrome plated brass tube binaurals

Description

Spare part for Laubry® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009526	HS Code	90189020009
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		