

NON CHILL RING (RED) FOR BELL OF DUAL PULSE II SPENGLER



Ref. 507945

Subtle thermal comfort for your Dual Pulse II exams, right down to the smallest detail.

Description

Spare part for SPENGLER DUAL PULSE® II stethoscopes.
Coral-colored anti-cold ring for small bell.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030094	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		