

NON CHILL RING (BLACK) FOR CARDIOPRESTIGE II STETHOSCOPE SPENGLER



Ref. 507972



29/10/2025 01:59:00 – For information purposes only, can be subject to change

NON CHILL RING (BLACK) FOR CARDIOPRESTIGE II STETHOSCOPE SPENGLER



Ref. 507972

A discreet and effective anti-chill ring designed to enhance examination comfort with the Cardio Prestige II.

Description

Spare part for SPENGLER Cardio Prestige® II stethoscopes.
Anti-cold ring for adult roof in black.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030315	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		