

FLOATING DIAPHRAGM (SMALL) FOR MAGISTER II DUAL SPENGLER

Ref. 507975S



A floating fiberglass diaphragm designed to enhance auscultation with the Magister II Dual.

Description

Spare part for stethoscopes MAGISTER® II from SPENGLER.
Multifrequency fiberglass membrane for pediatric pavilion.

29/10/2025 02:08:15 – For information purposes only, can be subject to change

FLOATING DIAPHRAGM (SMALL) FOR MAGISTER II DUAL SPENGLER



Ref. 507975S

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033187	HS Code	90181990
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.002

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---