

CHARGEUR ALIMENTATION EU ET CORDON POUR SPOT-CHECK 2VS



Ref. 546023

Reliable power supply to keep your Spot-Check 2VS running smoothly every day.

Description

European power charger and cord for SPENGLER spot-check station 2VS

Features

CE Marking	NONE	Medical device	No
Bar code	3700446029692	HS Code	85049099
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.093

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---