

INFRARED THERMOMETER FOR 4VS MONITOR



Ref. 551644

A reliable ear thermometer designed to fit seamlessly with the Spengler 4VS monitor.

Description

Infrared Ear Replacement Thermometer for Vital Signs Monitor SPENGLER 4VS

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033132	HS Code	9025192090
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.162

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---