

REPLACEMENT HANDLE FOR SMARTLED 5500 BLUEBERRY

Ref. 570519B



Intended for medical lighting use.

Description

A vibrant shade designed for better visibility and daily use with the Smartled 5500 otoscope.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446041748	HS Code	90189042
Country	Pakistan		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.038
Bare Unit / Weight (kg)	0.038		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		