

INSUFFLATION BULB FOR PNEUMATIC TEST

Ref. A5.201.112



29/10/2025 02:06:27 – For information purposes only, can be subject to change

INSUFFLATION BULB FOR PNEUMATIC TEST

Ref. A5.201.112

Enables diagnosis of movement/flexibility of the eardrum.

Compatible with Luxamed otoscopes or the Spengler X Luxamed edition otoscope.

Description

The Spengler inflation bulb is designed for quick and accurate manual cuff inflation. Its ergonomic grip ensures ease of use in both clinical and mobile settings.

Features

CE Marking	NONE	Medical device	Class I
Bar code	04260381960225	HS Code	9018908400
Country	Germany		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES