

BLUE NYLON ARMBAND, ADULT 31-40CM, 1T WITHOUT FITTING

SUNTECH

Ref. DSU230822

A durable and comfortable cuff designed for accurate blood pressure measurement in adults.

Description

Blue nylon adult cuff, designed for arm circumferences from 31 to 40 cm, compatible with direct-connect blood pressure monitors without adapters. Durable, comfortable, and easy to clean.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040219592	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 29/10/2025 01:49:16 – For information purposes only, can be subject to change