


MODULE TEMPERATURE COVIDIEN FILAC 3000 FOR STATION SUNTECH CT40

SUNTECH

Ref. DSU551403



CE
0413



Accurate and fast temperature readings with the Filac 3000 module, seamlessly compatible with the CT40 station.

Description

Covidien Fliac 3000 Tempereure Module for Suntech CT40 Station

29/10/2025 02:07:36 – For information purposes only, can be subject to change

MODULE TEMPERATURE COVIDIEN FILAC 3000 FOR STATION SUNTECH CT40

SUNTECH

Ref. DSU551403

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446032432	HS Code	9018.19.5500
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	23	Packed unit / Height (cm)	14
Packed unit / Length (cm)	19	Packed unit / Weight	0.612
Bare Unit / Weight (kg)	0.34		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---