

PRINTER PAPER CT50

Ref. [DSU222209](#)

SUNTECH

Thermal rolls designed to ensure clear and consistent printing with your SUNTECH CT50 monitor.

Description

Set of 4 thermal paper rolls compatible with the SUNTECH CT50 monitor, designed for clear and long-lasting printouts of medical records.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040214276	HS Code	9018199560
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.188

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---