

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE CHILD (12-19CM) - BAYONET CONNECTOR

SUNTECH

Ref. [DSU231011](#)



CE



A reusable pediatric cuff designed for comfortable and accurate readings with Bayonet connector.

Description

Reusable nylon cuff sold by 5. Single tube, bayonet coupling.

29/10/2025 02:02:11 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE CHILD (12-19CM) - BAYONET CONNECTOR

SUNTECH

Ref. DSU231011

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935104145	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.231

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---