

SET OF 20 ADULT 23-33CM

Ref. DSU231065

SUNTECH



CE



A single-patient adult cuff designed for reliable readings with Bayonet connector.

Description

Set of 20 Suntech "SOFT" single-patient armbands, royal blue.
Adult size 23-33cm, rectus bayonnette connection.

29/10/2025 02:01:34 – For information purposes only, can be subject to change

SET OF 20 ADULT 23-33CM

Ref. DSU231065

SUNTECH

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014179	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.319

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---