

**SET OF 20 SUNTECH MONOBLOC
ARMBANDS DOUBLE NEONATAL SINGLE
PATIENT TUBING N°4 (7-13CM)****SUNTECH****Ref. DSU231108**

Cuffs designed for the tiniest patients, ensuring reliable and gentle monitoring in neonatal care.

Description

This pack of 20 Suntech monobloc double-tube cuffs is tailored for neonatal use (Size 4 – 7 to 13 cm). It ensures accurate blood pressure readings while providing optimal comfort for the most delicate patients.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014872	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.433

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---