



SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

Ref. DSU250899

SUNTECH



CE
0413



A rigid sensor designed for accurate SpO₂ monitoring in adults, compatible with Suntech 247 and CT40 devices.

Description

SpO₂ articulated finger sensor for adult cable with a length of 3 meters for Suntech CT 40 station

29/10/2025 02:01:20 – For information purposes only, can be subject to change

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

SUNTECH

Ref. DSU250899

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	0108698870322340	HS Code	90181990
Warranty (months)	3	Country	Turkey
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	0.115
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		