


SPO2 REUSABLE DIGITAL SENSOR FOR ADULT

Ref. DSU250900

SUNTECH



CE
0413



Reliable and comfortable oxygen monitoring, designed for SunTech devices.

Description

SpO2 articulated finger sensor for adults with cable with a length of 3 meters for Suntech 247 station

29/10/2025 02:11:08 – For information purposes only, can be subject to change

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT

SUNTECH

Ref. DSU250900

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	010869887032265	HS Code	90181990
Country	Turkey		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.052
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		