

## REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

Ref. M1.121.513

**LUXAMED®**



CE



## REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

**LUXAMED®**Ref. M1.121.513

Reflex Hammer according to Buck | small

### Description

- Reflex hammer made of aluminium
- Optimal weight distribution
- Powder coating
- With neurological sensor tip
- Length: 18 cm | weight: 90 g
- Without brush due to hygienic reasons

Available in :

White: M1.121.914

Blue :M1.121.812

Pink :M1.121.612

Yellow: M1.121.513

Kid : M1.131.914

Black M1.121.112

Green: M1.121.412

Purple: M1.121.312

## REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

**LUXAMED®**

Ref. M1.121.513

### Features

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | CE             | Medical device | Class I  |
| Bar code          | 04260381960263 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Germany  |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

#### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature  | 40 |
| Minimum humidity    | 10  | Maximum humidity | 90 |

#### OTHER FEATURES

|        |        |
|--------|--------|
| Colors | YELLOW |
|--------|--------|