

SKY BLUE STAINLESS STEEL SINGLE TUBE BINAURALS

Ref. 507901



sky blue stainless steel single tube binaurals

Description

Spare part for Magister® stethoscopes

📅 29/10/2025 20:28:45 – For information purposes only, can be subject to change

SKY BLUE STAINLESS STEEL SINGLE TUBE BINAURALS

Ref. 507901



Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009434	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		