

PEARL MAGENTA CHROME PLATED BRASS SINGLE TUBE BINAURALS



Ref. 507915

pearl magenta chrome plated brass single tube binaurals

Description

Spare part for Pulse® & Dual Pulse® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446003890	HS Code	90181990
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---