

LED BULB FOR LITESTICK

Ref. 570690



Deliver clear and consistent lighting for medical exams with a reliable and long-lasting LED bulb.

Description

Spare part for Litestick® penlight.

📅 29/10/2025 20:25:35 – For information purposes only, can be subject to change

LED BULB FOR LITESTICK

Ref. 570690



Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009908	HS Code	90189084
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		