

ADU/PED/NEO NIBP TUBE (WHITE) 3M FOR NC3/NC5

COMEN

Ref. [BCBRAS12](#)

3-meter cuff extension – compatible with Comen NC3/NC5, flexible and extended reach.

Description

This 3-meter extension connects blood pressure cuffs to Comen NC3/NC5 monitors, providing greater mobility for the patient. Ideal for clinical settings where flexibility is key.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040219332	HS Code	9018390000
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---