

**FLOATING DIAPHRAGM/GREY RETAINING
RING ADULT SIZE CHESTPIECE****Ref. 507931**

Adult multi-frequency diaphragm – precise auscultation at varying pressures, with grey ring for Spengler stethoscope.

Description

Spare part for Cardio Prestige® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009625	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		