

NON CHILL RING (RED) FOR DUAL PULSE II SPENGLER



Ref. 507945B

Coral anti-chill ring – enhanced comfort for the large bell of the Dual Pulse II, ideal for extended examinations.

Description

Spare part for SPENGLER DUAL PULSE® II stethoscopes.
Coral color anti-cold ring for large pavilion.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030650	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.002
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		