

FLOATING DIAPHRAGM FOR CARDIOPRESTIGE II SPENGLER

Ref. 507976



Floating fiberglass diaphragm – enhanced acoustic precision for Cardio Prestige II stethoscope.

Description

Spare part for SPENGLER Cardio Prestige® II stethoscopes.
Multi-frequency fiberglass membrane for adult pavilion.

📅 29/10/2025 20:23:59 – For information purposes only, can be subject to change

FLOATING DIAPHRAGM FOR CARDIOPRESTIGE II SPENGLER



Ref. 507976

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030353	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		