

FLOATING DIAPHRAGM (SMALL) FOR CARDIOPRESTIGE II SPENGLER

Ref. 507977



FLOATING DIAPHRAGM (SMALL) FOR CARDIOPRESTIGE II SPENGLER



Ref. 507977

Floating fiberglass diaphragm – designed for the small bell of the Cardio Prestige, offering acoustic precision and listening comfort.

Description

Spare part for SPENGLER Cardio Prestige® II stethoscopes.
Multifrequency fiberglass membrane for pediatric pavilion.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030360	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		