

WHEEL STAND FOR SPOT-CHECK 2VS

Ref. 546024



Rolling stand – smooth mobility and enhanced stability for the Spengler 2VS spot-check station.

Description

Roller base for SPENGLER spot-check station 2VS

📅 29/10/2025 20:15:10 – For information purposes only, can be subject to change

WHEEL STAND FOR SPOT-CHECK 2VS

Ref. 546024



Features

CE Marking	NONE	Medical device	No
Bar code	3700446029715	HS Code	94029012003
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	10.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		