


SPO2 PEDIATRIC SENSOR FOR 4VS MONITOR

Ref. 551631



CE
0537



Pediatric SpO₂ sensor – accurate readings for young patients, designed for Spengler 4VS monitor.

Description

SPO2 Pediatric Vital Signs Monitor 4 VS SPENGLER

29/10/2025 20:20:13 – For information purposes only, can be subject to change

SPO2 PEDIATRIC SENSOR FOR 4VS MONITOR

Ref. 551631



Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031183	HS Code	90181990
Warranty (months)	6	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.093

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---