


SPO2 NEONATE SENSOR FOR 4VS MONITOR

Ref. 551632



CE
0537



Neonatal SpO₂ sensor – gentle and accurate monitoring from day one, designed for Spengler 4VS monitor.

Description

SPO2 Neonate Sensor for SPENGLER 4 VS Vital Signs Monitor

29/10/2025 20:26:38 – For information purposes only, can be subject to change

SPO2 NEONATE SENSOR FOR 4VS MONITOR

Ref. 551632



Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031190	HS Code	90181990
Warranty (months)	6	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.086

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---