

REPLACEMENT HANDLE FOR SMARTLED 5500 PEDIATRIC



Ref. 570519P

Intended for medical lighting use.

Description

Ergonomic replacement handle in mint green, designed for the Smartled 5500 pediatric otoscope. Provides a secure grip and easy identification in clinical settings.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446041724	HS Code	90189042
Country	Pakistan		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.038
Bare Unit / Weight (kg)	0.038		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--