

SPIROMETER DATOSPIR TOUCH EASY T

Ref. DSSPIR02

Sibelmed



CE
0197

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
1 YEAR
1
GUARANTEE

SPIROMETER DATOSPIR TOUCH EASY T

Ref. DSSPIR02



Assess respiratory function accurately with a user-friendly touch interface.

Description

The DATOSPIR TOUCH EASY T is a compact and intuitive electronic spirometer designed for primary care, occupational medicine, and pulmonary function units. It features a touchscreen interface, rechargeable battery, integrated weather station and automatic quality control to ensure reliable and accurate respiratory assessments.

Features

CE Marking	CE0197	Medical device	Class IIa
HS Code	90200090	Warranty (months)	12
Country	Spain		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		