


PLASTIC DISPOSABLE NOSE CLIP SIBELMED (X100)

Ref. DSSPIR03



CE



PLASTIC DISPOSABLE NOSE CLIP SIBELMED (X100)

Ref. [DSSPIR03](#)

Ensure reliable respiratory measurements with optimal nasal comfort and guaranteed hygiene.

Description

Designed for spirometry tests, these single-use nose clips prevent air leakage through the nostrils, ensuring accurate measurements. Their ergonomic design with soft pads and a flexible frame provides optimal patient comfort. Each clip is individually wrapped for maximum hygiene.

Features

CE Marking	CE	Medical device	Class I
Bar code	8699089770403	HS Code	90192090
Warranty (months)	12	Country	Turkey

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	375
Bare Unit / Weight (kg)	310		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES