

ADULT 23-33CM X20

Ref. DSU231085

SUNTECH



CE



Ensure reliable blood pressure measurement with a comfortable cuff designed for single-patient use.

Description

Set of 20 Suntech "SOFT" single-patient armbands, royal blue.
Double tubing, Adult size 23-33 cm, screw connection.

29/10/2025 20:22:12 – For information purposes only, can be subject to change

ADULT 23-33CM X20

Ref. DSU231085

SUNTECH

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014520	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	16
Packed unit / Length (cm)	25	Packed unit / Weight	1.444
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		