

MASIMO SET SPO2 MODULE SUNTECH CT40

SUNTECH

Ref. [DSU551402](#)

Masimo SET SpO₂ module – advanced technology for reliable and continuous monitoring on Suntech CT40 station.

Description

Extra-large single-use cuffs designed for adult patients requiring extended sizing. The Bayonet connector ensures quick and secure attachment to compatible devices.

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3661040217700	HS Code	9018.19.5500
Warranty (months)	24		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--