


ELITE BAGS BRIEFCASE - COMMUNITY - BLUE

Ref. MECOMM02_RC

 **ELITE BAGS**



 30/10/2025 21:03:07 – For information purposes only, can be subject to change

ELITE BAGS BRIEFCASE - COMMUNITY - BLUE



Ref. MECOMM02_RC

Description

It is the ideal compact case for professionals who provide limited medical assistance, such as nurses, or carry very specific equipment. With an interior design that allows to modify the configuration of the spaces in order to accommodate different types of materials.

- A large front pocket for personal items (gloves, needle collector ...)
- Space with extractable separators
- Isothermal amp maker with 19 elastics
- Removable shoulder straps

- Padlock closure
- Strip for positioning the trolley (optional)- Measures: 35 x 25.5 x 14 cm
- Capacity: 12.5 L
- Weight: 1.35 kg

- Material: 420D Serge polyester jacquard. Washable
- Maximum recommended capacity: 3kg

- Sold empty
- Available in blue (MECOMM02_B) and pink (MECOMM02_RC)

ELITE BAGS BRIEFCASE - COMMUNITY - BLUE



Ref. [MECOMM02_RC](#)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564541789	HS Code	42029270
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	25.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	35
Packed unit / Weight	1.5	Bare Unit / Weight (kg)	1.35

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Pink
--------	------