

EMERGENCY BACKPACK - TRAPEZOIDAL - BLACK (EM13.019)

Ref. METRAP01_R

ELITE BAGS



30/10/2025 20:48:27 – For information purposes only, can be subject to change

EMERGENCY BACKPACK - TRAPEZOIDAL - BLACK (EM13.019)



Ref. METRAP01_R

Trapezoidal emergency bag, designed for rescues and interventions in different scenarios, such as in difficult access areas

Description

- Plastic surface (waterproof and easy to clean)
- 4 removable compartments with transparent upper part. With handles of different colors (red, blue, green, yellow) to differentiate the content.
- 3 net pockets
- Outside, 3 pockets for first aid equipment (gloves, bandages, etc.)
- Reflective strips.
- U-shaped zipper, with double slider for quick opening
- 4 plastic loops for more security
- Top and side handle
- Reinforced and protected base
- Hidden backpack straps
- Outdoor measures: 55 x 28 x 20 cm / 30.8 L
- Weight: 1.62 kg
- Material: 600D Polyester with a PU coating
- Color: Black or Red
- Max. recommended capacity: 5 kg
- Sold empty

EMERGENCY BACKPACK - TRAPEZOIDAL - BLACK (EM13.019)



Ref. [METRAP01_R](#)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542038	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	39	Bare Unit / Width (cm)	28
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	53	Bare Unit / Length (cm)	55
Packed unit / Weight	1.64	Bare Unit / Weight (kg)	1.62

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Red	Packed unit / Quantity	1
--------	-----	------------------------	---