


SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC FINGER CLIP FOR MP3



Ref. 222261



CE
0537



31/10/2025 21:54:56 – For information purposes only, can be subject to change

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC FINGER CLIP FOR MP3



Ref. 222261

Reliable SpO₂ monitoring for children with OxiMax technology.

Description

NELCORR OXIMAX compatible spo2 sensor, pediatric for oximeter MASTER PALM 3

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030148	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.094

SPECIFIC STORAGE CONDITIONS

Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---