

EMERALD STAINLESS STEEL SINGLE BINAURALS

Ref. 507906



emerald stainless steel single binaurals

Description

Spare part for Magister® stethoscopes

31/10/2025 21:30:15 – For information purposes only, can be subject to change

EMERALD STAINLESS STEEL SINGLE BINAURALS

Ref. 507906



Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009472	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.065
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		