

# YELLOW CHROME PLATED BRASS SINGLE TUBE BINAURALS



Ref. 507919

yellow chrome plated brass single tube binaurals

## Description

Spare part for Laubry® stethoscopes

## Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | NONE          | Medical device            | Class I  |
| Bar code                           | 3700446009564 | HS Code                   | 90181990 |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.08     |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | NA            | Max temperature           | NA       |
| Minimum humidity                   | NA            | Maximum humidity          | NA       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Packed unit / Quantity             | 1             |                           |          |