

FLOATING DIAPHRAGM/BLACK RETAINING RING FOR CHILDREN SIZE CHESTPIECE

Ref. 507932



Precise auscultation of high and low frequencies tailored for pediatric assessments.

Description

Spare part for Cardio Prestige® stethoscopes

31/10/2025 21:41:42 – For information purposes only, can be subject to change

FLOATING DIAPHRAGM/BLACK RETAINING RING FOR CHILDREN SIZE CHESTPIECE



Ref. 507932

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009632	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		