

FLOATING DIAPHRAGM/GREY RETAINING RING FOR CHILDREN SIZE CHESTPIECE

Ref. 507933



Fine auscultation of medical sounds with a diaphragm tailored for pediatric use.

Description

Spare part for Cardio Prestige® stethoscopes

31/10/2025 21:45:09 – For information purposes only, can be subject to change

FLOATING DIAPHRAGM/GREY RETAINING RING FOR CHILDREN SIZE CHESTPIECE



Ref. 507933

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009649	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		