

FLOATING DIAPHRAGM/BLACK RETAINING RING CARDIO SONIC



Ref. 507934

A multi-frequency diaphragm for precise and comfortable auscultation with the Cardio Sonic stethoscope.

Description

Spare part for Cardio Sonic® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009656	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		