

GREY RING FOR MAGISTER/PULSE BELL

Ref. 507940



A grey antifreeze ring for optimal comfort with the bell chestpiece of the Magister Dual Pulse stethoscope.

Description

Spare part for Pulse® & Dual Pulse® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009694	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		