

NON CHILL RING (PINK) FOR BELL OF DUAL PULSE II SPENGLER



Ref. 507947

Thermal comfort for examinations with the small bell of the Dual Pulse II stethoscope.

Description

Spare part for SPENGLER DUAL PULSE® II stethoscopes.
Anti-cold ring black pink small bell.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030117	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		