


SPO2 SENSOR PEDIATRIC FOR SPOT-CHECK 2VS

Ref. 546011



CE
0123



Accurate and child-friendly monitoring for reliable readings with the Spot-check 2VS.

Description

SPO2 pediatric sensor for SPENGLER spot-check station 2VS.

31/10/2025 21:50:56 – For information purposes only, can be subject to change

SPO2 SENSOR PEDIATRIC FOR SPOT-CHECK 2VS

Ref. 546011



Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446029616	HS Code	90181990
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.056

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---