


## SPO2 SENSOR NEONATE FOR SPOT-CHECK 2VS

Ref. 546012



CE  
0123



Gentle and secure SpO<sub>2</sub> monitoring for newborns with the Spot-check 2VS.

### Description

SPO2 neonate sensor for SPENGLER spot-check station 2VS.

31/10/2025 21:54:38 – For information purposes only, can be subject to change

## SPO2 SENSOR NEONATE FOR SPOT-CHECK 2VS

Ref. 546012



### Features

|                   |               |                |           |
|-------------------|---------------|----------------|-----------|
| CE Marking        | CE0123        | Medical device | Class IIb |
| Bar code          | 3700446029623 | HS Code        | 90181990  |
| Warranty (months) | 6             | Country        | China     |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |      |
|---------------------------|----|---------------------------|------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc   |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.07 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|