

NEONATE CUFF 6-11CM FOR SPOT-CHECK 2VS

Ref. 546016



A cuff designed for newborns, ensuring accurate readings with the Spot-Check 2VS Spengler.

Description

Neonate cuffs (6-11cm) for spot-check station 2VS from SPENGLER.

31/10/2025 21:52:26 – For information purposes only, can be subject to change

NEONATE CUFF 6-11CM FOR SPOT-CHECK 2VS

Ref. 546016



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029630	HS Code	9018901000
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.029

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---