

PEDIATRIC CUFF 10-19CM FOR SPOT-CHECK 2VS

Ref. 546017



Tailored for young patients with fast and accurate spot-check measurement.

Description

Pediatric armbands (10-19cm) for spot-check station 2VS from SPENGLER.

31/10/2025 21:40:29 – For information purposes only, can be subject to change

PEDIATRIC CUFF 10-19CM FOR SPOT-CHECK 2VS

Ref. 546017



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029647	HS Code	9018901000
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.037

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---