

WALL BRACKET FOR 4VS MONITOR

Ref. 551611



A sturdy and discreet wall mount to keep the Spengler 4VS monitor within easy reach.

Description

This wall mount is designed to securely hold the Spengler 4VS monitor. It offers a fixed installation that optimizes workspace and ensures easy access during medical assessments.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040213835	HS Code	85049099
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.58
Bare Unit / Weight (kg)	3.58		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

31/10/2025 21:47:46 – For information purposes only, can be subject to change